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EN | INSTRUCTIONS FOR USE

DESCRIPTION

Kysense® Precise / Kysense® Touch / Kysense® Define / Kysense® Volume / Kysense® Extreme is a sterile, transparent and resorbing gel of cross-linked hyaluronic acid of non-animal origin. The product is supplied in a plastic syringe and it is intended for single use only.

The volume of each syringe is 1 mL for **Kysense® Precise / Kysense® Define / Kysense® Volume / Kysense® Extreme**. **Kysense® Touch** exists with a volume of 0.5 mL per syringe.

Each box contains four graduated pre-filled syringes of **Kysense® Precise** or **Kysense® Touch** or **Kysense® Define** or **Kysense® Volume** or **Kysense® Extreme**, an instruction leaflet, a set of traceability labels and four single-use sterile needles of 30G ½ for **Kysense® Precise / Kysense® Touch / Kysense® Define** and of 27G ½ for **Kysense® Volume / Kysense® Extreme**. For each syringe, one label is to be given to the patient, the other one is to be added to the patient record.

The content of **Kysense® Precise / Kysense® Touch / Kysense® Define / Kysense® Volume / Kysense® Extreme** syringes is sterilised by moist heat.

The 30G ½ or 27G ½ needles are sterilised by radiation.

The intended purpose of **Kysense® Precise / Kysense® Touch / Kysense® Define / Kysense® Volume / Kysense® Extreme** is to modify the anatomy and / or to alleviate a physiological process at the level of the face.

COMPOSITION

Kysense® Touch Kysense® Precise	<ul style="list-style-type: none"> • Cross-linked* hyaluronic acid: • Phosphate buffer saline One syringe contains 1 mL of Kysense® Precise or 0.5 mL of Kysense® Touch	15 mg / mL q.s. 1 mL
Kysense® Define	<ul style="list-style-type: none"> • Cross-linked* hyaluronic acid: • Phosphate buffer saline One syringe contains 1 mL of Kysense® Define	18 mg / mL q.s. 1 mL
Kysense® Volume	<ul style="list-style-type: none"> • Cross-linked* hyaluronic acid: • Phosphate buffer saline One syringe contains 1 mL of Kysense® Volume	21 mg / mL q.s. 1 mL
Kysense® Extreme	<ul style="list-style-type: none"> • Cross-linked* hyaluronic acid: • Phosphate buffer saline One syringe contains 1 mL of Kysense® Extreme	24 mg / mL q.s. 1 mL

*The hyaluronic acid has been crosslinked with BDDE.

INTENDED USES

Kysense® Touch Kysense® Precise	<p>Kysense® Precise / Touch is an injectable gel intended for:</p> <ul style="list-style-type: none"> • Facial reconstruction of structural defects from congenital or medical origin: scar tissue treatment • Aesthetic treatment of fine lines, medium-sized skin depressions and lips definition or enhancement at the level of the face. <p>Kysense® Precise / Touch is intended to be administrated via dermis injection or lips mucosa injection by an authorized health professional.</p>
Kysense® Define	<p>Kysense® Define is an injectable gel intended for:</p> <ul style="list-style-type: none"> • Facial reconstruction of structural defects from congenital or medical origin: scar tissue treatment • Aesthetic treatment of deep skin depressions and lips enhancement at the level of the face. <p>Kysense® Define is intended to be administrated via deep dermis injection or lips mucosa injection by an authorized health professional.</p>
Kysense® Volume	<p>Kysense® Volume is an injectable gel intended for:</p> <ul style="list-style-type: none"> • Facial reconstruction of structural defects from congenital or medical origin: volume lost by HIV-associated lipoatrophy • Aesthetic treatment of deep skin depressions at the level of the face and augmentation of the volume of facial tissues. <p>Kysense® Volume is intended to be administrated via injection into the subcutaneous fat tissue or into the suprapariosteal zone by an authorized health professional.</p>
Kysense® Extreme	<p>Kysense® Extreme is an injectable gel intended for:</p> <ul style="list-style-type: none"> • Facial reconstruction of structural defects from congenital or medical origin: volume lost by HIV-associated lipoatrophy • Aesthetic treatment of facial volume restoration. <p>Kysense® Extreme is intended to be administrated via injection into the subcutaneous fat tissue or into the suprapariosteal zone by an authorized health professional.</p>

MODE OF ACTION

Kysense® Precise / Kysense® Touch / Kysense® Define / Kysense® Volume / Kysense® Extreme creates a volume in the tissue that corrects facial anatomy.

Kysense® Precise / Kysense® Touch / Kysense® Define / Kysense® Volume / Kysense® Extreme is intended to be wholly resorbed over time, which progressively reduce the volumizer effect of the product at the injection site. The kinetic of this absorption depends on several factors among which the injected quantity, the depth of injection and the patient metabolism.

CONTRA-INDICATIONS

- Do not inject patients having known risk of hypersensitivity to hyaluronic acid.
- Do not inject patient having history of severe allergy or anaphylactic shock.
- Do not inject patients suffering from autoimmune disease.
- Do not inject patients having cutaneous disorders, inflammation or infection (herpes, acne, etc.) at the treatment site or nearby.
- Do not inject patients for whom the medical history shows a sensitivity that could lead to a reaction to the treatment.
- Do not use in patients with bleeding disorders or in patients who are undergoing treatment with thrombolytics or anticoagulants.
- Do not inject patients with a tendency to form keloids, hypertrophic scars or any other healing disorders.
- Do not use during pregnancy or breast-feeding.
- Do not use for children.

WARNING

- Do not inject into blood vessels (intra-vascular) or directly near blood vessels (for avoiding vessels compression).
- Do not inject in area other than the face (e.g. hands, body).
- Do not inject **Kysense® Precise / Kysense® Touch** into the eyelids. The application of **Kysense® Precise / Kysense® Touch** in the periorbital area is reserved to specialists specifically trained in this technique and having sound knowledge of the physiology for this particular area.
- Do not inject **Kysense® Define / Kysense® Volume / Kysense® Extreme** in the periorbital area (eyelids, palpebro-malar groove, crow's feet) and glabellar region.
- Intra-muscular injection is not recommended.
- There is no available clinical data about injection of **Kysense® Precise / Kysense® Touch / Kysense® Define / Kysense® Volume / Kysense® Extreme** into an area which has already been treated with another brand of resorbing product. It is therefore recommended not to inject in areas in which another brand of resorbing product is not still fully resorbed.
- Do not inject jointly with another brand of resorbing product in the same area.
- Do not inject jointly with a permanent implant.
- Do not mix with other products.
- Do not re-sterilize the product (the syringe and / or the needle).
- Do not inject more than 20 mL per year per patient.
- Do not inject more than 2 mL per treatment area during each session.
- Do not overcorrect.
- Do not use after expiry date or if the packaging is damaged or the syringe opened.
- Do not re-use (risk of contamination). The product is intended for single use only.
- Do not use a syringe for several patients (risk of cross-contamination).
- Do not use simultaneously with chemical peels, dermabrasion, laser treatment or radiofrequency treatment.

PRECAUTIONS FOR USE

- Patients must be questioned regarding their medical history and must be informed about expected outcome, warning, contra-indications, precautions and potential side effects.
- The quantity to inject depends on the area to treat and the expectation of the patient. Patient with unrealistic expectation must not be treated.
- Aseptic techniques and standards practices are to be observed.
- The content of the syringe is sterile until opening. External surface of the syringe is not sterile and should not be used in a controlled environment such as operating room. If needed, the surface of the syringe must be cleaned up in respect of the aseptic techniques in place.
- The patients on medications that affect platelet functions (aspirin or non-inflammatory drugs) must be warned of the potential increased risks of haematomas and bleeding during injection.
- Special caution must be exercised when treating areas in close proximity to vulnerable structures such as vessels and nerves.
- Advise the patients to avoid any extended exposure to the sun, UV rays and temperatures below 0°C, as well as sauna or hammam sessions during the two weeks following the injection procedure.
- Recommend the patient not use cosmetic products (e.g. makeup, skin care products) during the 12 hours following the injection treatment.
- Recommend the patient to avoid applying intense pressure or massaging the treatment site for a few days following the injection.

ASSEMBLY OF NEEDLE TO SYRINGE

- Use of the supplied 30G ½ (for **Kysense® Precise / Kysense® Touch / Kysense® Define**) and 27G ½ (for **Kysense® Volume / Kysense® Extreme**) sterile needle is recommended. Use of another needle could not ensure the correct fitting to the syringe and / or the correct delivery of the product.
- Remove the tip cap of the syringe.
- Grasp the hub.
- Press and turn the needle clockwise on the syringe until it is tight.
- Make sure that the hub of the needle is as tightly fixed to the bottom of the luer-lock as possible and remove the shield of the needle just before injection.
- Failure to comply with these precautions could cause disengagement of the needle and / or product leakage at luer-lock level.

INTERACTIONS

Treatment with **Kysense® Precise** / **Kysense® Touch** / **Kysense® Define** / **Kysense® Volume** / **Kysense® Extreme** in combination with drugs and other devices has not been clinically evaluated.

The treatment with **Kysense® Precise** / **Kysense® Touch** / **Kysense® Define** / **Kysense® Volume** / **Kysense® Extreme** can be performed in combination with an appropriate quantity of lidocaine solution.

PROCEDURE OF ADMINISTRATION

Treatment with **Kysense® Precise** / **Kysense® Touch** / **Kysense® Define** / **Kysense® Volume** / **Kysense® Extreme** is only intended to be administered by an authorized healthcare professional in accordance with the local legislation.

The healthcare professional shall have a deep knowledge of the anatomy of the area to be treated and shall be used to performing similar intervention (i.e. cross-linked hyaluronic injection in the face area).

Kysense® Precise / **Kysense® Touch** / **Kysense® Define** / **Kysense® Volume** / **Kysense® Extreme** is to be used as supplied. Modification or use of the product outside of the instruction for use may adversely impact the sterility and / or the performance of the device.

A correct injection technique is important for the clinical outcome of the treatment. The injection technique varies depending on the treatment site, the depth of injection and the gel quantity to be injected.

- 1) Before starting treatment, the patient should be informed of the product's indications, contra-indications, precautions, warnings, expected results and potential adverse events and questioned on her / his medical history.
- 2) The injected area must be disinfected with an appropriate antiseptic solution.
- 3) Assemble the supplied sterile needle and syringe.
- 4) Do not apply extensive pressure to the syringe at any time. If such resistance is encountered, remove and replace the needle and check for function.
- 5) **Kysense® Precise** / **Kysense® Touch** / **Kysense® Define** / **Kysense® Volume** / **Kysense® Extreme** should be administered slowly in the tissues. The injection volume varies depending on the correction required. Do not overcorrect.
Kysense® Precise / **Kysense® Touch** is indicated for injection into the dermis or lips mucosa.
Kysense® Define is indicated for injection into the deep dermis or lips mucosa.
Kysense® Volume is indicated for injection into the subcutaneous fat tissue or into the supraperiosteic zone.
Kysense® Extreme is indicated for injection into the subcutaneous fat tissue or into the supraperiosteic zone.
- 6) For optimal clinical results, gentle massage by the healthcare professional of the treated area is recommended to ensure uniform **distribution of the material**.

ADVERSE EVENTS

Healthcare professional shall inform the patients about potential adverse events related to the injection procedure and / or the use of the product **Kysense® Precise** / **Kysense® Touch** / **Kysense® Define** / **Kysense® Volume** / **Kysense® Extreme**, occurring immediately after injection or after a delay.

These adverse events include but are not limited to:

-Events which are naturally resolved within one week in most cases:

- Injection-related events and / or inflammatory reactions such as bleeding, ecchymosis, erythema, haematoma, skin redness, bruising, swelling, oedema and infection which may be associated with local pain or itching, occurring after injection.
- Sensitivity at the injection site.
- Hardness, lump or nodule at the injection site.
- Skin coloration or discoloration at the injection site.

-Events which are delayed after the injection:

- Immediate or delayed hypersensitivity to hyaluronic acid.
- Nodule, abscesses or granuloma.
- Vascular compromise may occur due to inadvertent intravascular injection or as a result of vascular compression associated with implantation of any injectable soft tissue filler. This may manifest as blanching, discoloration, necrosis or ulceration at the implant site or in the area supplied by the blood vessels affected; or rarely as ischemic events in other organs due to embolization. Isolated rare cases of ischemic events affecting the eye leading to visual loss, and the brain resulting in cerebral infarction, following facial aesthetic treatments have been reported. These rare cases of vascular events are mostly reported in glabella, in and around the nose, in forehead and in periorbital area.
- Infection or reactivation of a previous infection.
- Displacement of the gel.
- Inflammatory reactions which persist for more than one week, or any other adverse event which develops, must be reported by the patients to their healthcare professional. In this case, if required, the healthcare professional should use an appropriate treatment.

Any adverse events and complications must be reported to the local representative of Kylane Laboratoires.


WASTE DISPOSAL

- At the end of the session, opened and used syringe as well as all injection materials (e.g. needle) shall be discarded in respect of the medical practices and in accordance with the applicable regulations on waste disposal (e.g. specific waste container).

STORAGE CONDITIONS

- Shelf life of the product is indicated on each syringe and product must not be used after this date.
- Store between 2°C and 25°C away from direct sunlight.
- Make sure there is no visible signs of damage to the packaging before use.

MEANING OF THE SYMBOLS ON THE PACKAGING:

LOGO	SIGNIFICATION	LOGO	SIGNIFICATION
	Batch code		Sterile fluid path - Moist heat (the content of the syringe (gel) is sterile)
	Use-by date		Do not use if package is damaged
	Sterilised using moist heat (syringes)		Keep away from sunlight
	Do not re-use		Fragile; handle with care
	Caution		Temperature limit (+2°C - +25°C)
	Manufacturer		Read instructions for use
	CE-mark 2797 refers to the Notified Body		Sterilised using irradiation (needles)

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Distributed by: Circa Skin Ltd.
Sinclair France SAS